

Technical Datasheet

		EN ISO 13485 certified
Article - No:	02-062-0310	Revision: 4
Remarks		single use only

Specification

Product	Micro tubes PP, 1,1ml, Ø8,8x45 mm, 12x8-strips/rack, transparent, np pcr ready, autoclavable up to 121°C
Dimensions	length: 45,00 mm (+/- 0,2 mm) inner diameter: 8,8 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	1,1ml
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 75°C
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62 Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors free of human DNA free of RNase and DNase free of pyrogens
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 960 pcs per carton: 4.800 pcs per pallet: 144.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	