

Technical Datasheet

EN ISO 13485 certified		
Article - No:	02-062-1150	Revision: 2
Remarks		single use only

Specification

Product	Caps PE for micro tubes, strip with 8 caps, transparent, np pcr ready
Colour	transparent
Material Characteristics	polyethylene (PE)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 120
	pcs per carton: 600
	pcs per pallet: 36.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	