

Technical Datasheet

Date: 22.06.20

Page 1 of 1

EN ISO 13485 certified	
Article - No:	02-062-1800
Remarks	Revision: 2 single use only

Specification

Product	Micro tubes PP, 1,1ml, Ø8,8x45 mm, 8-strips, transparent, np pcr ready, autoclavable up to 121°C
Dimensions	length: 45,00 mm (+/- 0,2 mm)
	inner diameter: 8,8 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	1,1ml
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 75°C
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 120
	pcs per carton: 600
	pcs per pallet: 36.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	