

Technical Datasheet

EN ISO 13485 certified		
Article - No:	02-112-0001	Revision: 2
Remarks	Fitting plug: 02-633-0000	single use only

Specification

Product	Test tube PP, round bottom, 4ml, Ø12x55 mm, transparent, max. RCF 3.000g, autoclavable up to 121°C
Dimensions	outer diameter: 12,00 mm (+/- 0,2 mm)
	height: 55,00 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	4ml
Material Characteristics	polypropylene (PP)
Temperature range	-10°C - 121°C
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 3.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 1.000
	pcs per carton: 5.000
	pcs per pallet: 120.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

