

Technical Datasheet

EN ISO 13485 certified		
Article - No:	02-161-0001	Revision: 2
Remarks	Fitting plug: 02-633-0000	single use only

Specification

Product	Test tube PS, round bottom, 5ml, Ø13x75 mm, transparent, max. RCF 3.000g
Dimensions	Inner diameter: 12,73 mm (+/- 0,2 mm)
	Outer diameter: 13,00 mm (+/- 0,2 mm)
	Height : 74,55 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	5ml
Material Characteristics	polystyrene (PS)
Temperature range	-10°C – 70°C
Centrifugation	max. RCF 3.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 1.000
	pcs per carton: 4.000
	pcs per pallet: 80.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

