

Technical Datasheet

EN ISO 13485 certified		
Article - No:	02-221-0000	Revision: 1
Remarks	Fitting plug: 02-673-0000	single use only

Specification

Product	Test tube PS, round bottom, 12ml, Ø16x100 mm, transparent, max. RCF 3.000g
Dimensions	outer diameter: 16,10 mm (+/- 0,2 mm) height: 99,90 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	12ml
Material Characteristics	polystyrene (PS)
Temperature range	-10°C – 70°C
Centrifugation	max. RCF 3.000g
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 1.000 pcs per carton: 2.000 pcs per pallet: 40.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

nerbe plus GmbH & Co. KG
Porschestr. 25
D-21423 Winsen/Luhe
Tel: +49 (0)4171 7849 0
Fax: +49 (0)4171 6 45 45

info@nerbe-plus.de
www.nerbe-plus.de
HRA Lüneburg 203240
USt.-ID: DE 32583 0652

Pers. haftende Gesellschafterin:
NERBE R&D GmbH, Winsen/Luhe
GF: Philipp J.F. Nerbe
HRB Lüneburg 208826

Deutsche Bank AG
Kto. 0204 41200
IBAN: DE41 2407 0075 0020
4412 00
SWIFT / BIC DEUTDE2240

Hamburger Sparkasse
Kto. 1387 140831
IBAN: DE56 2005 0550 1387 1408 31
SWIFT / BIC HASPDE33XXX

