

Technical Datasheet

EN ISO 13485 certified		
Article - No:	02-252-0000	Revision: 3
Remarks	Fitting plug: 02-673-0000	single use only

Specification

Product	Test tube PP, round bottom, 20ml, Ø16x150 mm, transparent, max. RCF 3.000g, autocl. up to 121°C
Dimensions	outer diameter: 16,00 mm (+/- 0,2 mm)
	height: 150,00 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	20ml
Material Characteristics	polypropylene (PP)
Temperature range	-10°C - 121°C
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 3.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 850
	pcs per carton: 850
	pcs per pallet: 20.400
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

