

Technical Datasheet

		EN ISO 13485 certified
Article - No:	02-352-1000	Revision: 3
Remarks		single use only

Specification

Product	Test tube PP with enclosed screw cap PE, 12ml, Ø16x100 mm, round bottom, tube transparent, cap natural, grad., sterile, max. RCF 5.400g, autocl. up to 110°C
Dimensions	outer diameter: 16,00 mm (+/- 0,2 mm) height (without screw cap): 100,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap natural
Graduation	grad.
Nominal Capacity	12ml
Material Characteristics	tube: polypropylene (PP) cap: polyethylene (PE)
Temperature range	tube: -196°C up to +121°C cap: -50°C up to +80°C
Autoclavability	Autocl. up to 121°C (no cap)
Centrifugation	max. RCF 4.000g
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62 Food contact: raw material in compliance with regulation EC 1935/2004
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 1.000 pcs per carton: 1.000 pcs per pallet: 24.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	

