

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>02-381-8107</b>	<b>Revision: 3</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Test tube PS with fitted screw cap PE, 19ml, Ø16x150 mm, round bottom, tube transparent, cap red, sterile R
<b>Dimensions</b>	outer diameter: 16,00 mm (+/- 0,2 mm)
	height (without screw cap): 150,00 mm (+/- 0,2 mm)
<b>Colour</b>	tube transparent, cap red
<b>Nominal Capacity</b>	19ml
<b>Material Characteristics</b>	tube: polystyrene (PS)
	cap: polyethylene (PE)
<b>Temperature range</b>	-10°C – 70°C
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
<b>Purity Grade</b>	do not contain any pyrogens
<b>Sterilisation Process</b>	radiation sterilization SAL 10-6
<b>Shelf - Life</b>	36 months at controlled storage conditions
<b>Packaging</b>	pcs per bag: 100
	pcs per carton: 800
	pcs per pallet: 19.200
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
<b>Approved by nerbe plus quality department</b>	

