

# Technical Datasheet

EN ISO 13485 certified		
<b>Article - No:</b>	<b>02-401-0000</b>	<b>Revision: 3</b>
<b>Remarks</b>	<b>Fitting plug: 02-633-0000</b>	<b>single use only</b>

## Specification

<b>Product</b>	Test tube PS, conical bottom, 4ml, Ø11x75 mm, transparent, max. RCF 3.000g
<b>Dimensions</b>	outer diameter: 11,00 mm (+/- 0,2 mm)
	height: 75,00 mm (+/- 0,2 mm)
<b>Colour</b>	transparent
<b>Nominal Capacity</b>	4ml
<b>Material Characteristics</b>	polystyrene (PS)
<b>Temperature range</b>	-10°C – 70°C
<b>Centrifugation</b>	max. RCF 3.000g
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
<b>Packaging</b>	pcs per bag: 100
	pcs per carton: 4.000
	pcs per pallet: 80.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

**Approved by nerbe plus quality department**

00-000-0001/06.2020

