

# Technical Datasheet

|                        |  |
|------------------------|--|
| EN ISO 13485 certified |  |
| <b>Article - No:</b>   | <b>02-531-6007</b>                           |
| <b>Remarks</b>         | <b>Revision: 5</b><br><b>single use only</b> |

## Specification

|                                 |  |
|---------------------------------|--|
| <b>Product</b>                  | Test tube PS with fitted screw cap PE, 25ml, Ø25x90 mm, conical bottom, skirted, tube transparent, cap white, sterile R, max. RCF 4.000g |
| <b>Dimensions</b>               | tube:  |
|                                 | outer diameter: 25,00 mm (+/- 0,2 mm)  |
|                                 | height: 90,00 mm (+/- 0,2 mm)  |
|                                 | cap:   |
|                                 | outer diameter: 31,00 mm (+/- 0,2 mm)  |
|                                 | height: 17,00 mm (+/- 0,2 mm)  |
| <b>Colour</b>                   | tube transparent, cap white  |
| <b>Nominal Capacity</b>         | 25ml   |
| <b>Material Characteristics</b> | tube: polystyrene (PS)   |
|                                 | cap: polyethylene (PE)   |
| <b>Temperature range</b>        | -10°C – 70°C   |
| <b>Centrifugation</b>           | max. RCF 4.000g  |
| <b>Product Characteristics</b>  | biologically inert   |
|                                 | physiologically harmless   |
|                                 | no heavy metals according to EC/94/62  |
| <b>Purity Grade</b>             | do not contain any pyrogens  |
| <b>Sterilisation Process</b>    | radiation sterilization SAL 10-6   |
| <b>Shelf - Life</b>             | 60 months at controlled storage conditions   |
| <b>Packaging</b>                | pcs per bag: 100   |
|                                 | pcs per carton: 500  |
|                                 | pcs per pallet: 12.000   |
| <b>Storage</b>                  | dry condition: 40%-60% rel. humidity   |
|                                 | bags / cases orig. closed  |
|                                 | storage temp. +10°C to max. +30°C  |
|                                 | pal. not stackable   |



# Technical Datasheet

Date: 02.06.22

Page 2 of 2

Approved by nerbe plus quality department