

Technical Datasheet

		EN ISO 13485 certified
Article - No:	02-633-0000	Revision: 2
Remarks		single use only

Specification

Product	Ribbed plug PE, with high grip, with 2 sealing ribs, Ø 12 mm, natural
Dimensions	diameter: 12,00 mm (+/- 0,2mm)
Colour	natural
Material Characteristics	polyethylene (PE)
Temperature range	-50°C - 80°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 1.000
	pcs per carton: 20.000
	pcs per pallet: 480.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	