

Technical Datasheet

Date: 22.06.20

Page 1 of 1

		EN ISO 13485 certified
Article - No:	02-673-0000	Revision: 1
Remarks		single use only

Specification

Product	Ribbed plug PE, with high grip, with 2 sealing ribs, Ø16 mm, natural
Dimensions	diameter: 16,00 mm (+/- 0,2 mm)
	height: 21,00 mm (+/- 0,2 mm)
	Sealing ribs diameter: 14,50 mm (+/- 0,2 mm)
Colour	natural
Material Characteristics	polyethylene (PE)
Temperature range	-50°C - 80°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 1.000
	pcs per carton: 10.000
	pcs per pallet: 240.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	