

Technical Datasheet

		EN ISO 13485 certified
Article - No:	02-801-2000	Revision: 1
Remarks		single use only

Specification

Product	Cell counter vial PS, dust-free capped, 15ml, Ø22x45 mm, suitable for Labora, Contraves Hellige, Brand, tube transparent, cap natural
Dimensions	diameter: 22,00 mm (+/- 0,2 mm)
	height: 45,00 mm (+/- 0,2 mm)
Colour	tube transparent, cap natural
Nominal Capacity	15ml
Material Characteristics	polystyrene (PS)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Packaging	pcs per bag: 500
	pcs per carton: 2.000
	pcs per pallet: 24.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

