

Technical Datasheet

EN ISO 13485 certified		
Article - No:	02-908-4000	Revision: 1
Remarks	DMSO resistant, available in English on request	single use only

Specification

Product	Barcode label 50,8x5,6 mm, without print, for microtiter plates, white
Dimensions	length: 50,80 mm (+/- 0,2 mm)
	width: 5,60 mm (+/- 0,2 mm)
Colour	white
Material Characteristics	paper (PAP)
Packaging	pcs per carton: 100.000
	pcs per Rolle: 10.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	