

Technical Datasheet

EN ISO 13485 certified	
Article - No:	04-032-0350
Revision:	4
Remarks	single use only

Specification

Product	PCR microcentrifuge tube PP, 0,1ml, 8-strips, individually attached flat & transparent caps, low profile, highly transparent, np pcr ready, max. RCF 20.000g, autocl. up to 121°C
Colour	highly transparent
Nominal Capacity	0,1ml
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 121°C
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 20.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per inner cs: 120
	pcs per carton: 1.200
	pcs per pallet: 57.600
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

