

Technical Datasheet

EN ISO 13485 certified	
Article - No:	04-032-0800
Revision:	4
Remarks	single use only

Specification

Product	PCR microcentrifuge tube PP, 0,2ml, 8-strips, with attached stripe of 8 domed caps, transparent np pcr ready, max. RCF 6.000g, autocl. up to 121°C CE/IVD
Classification	CE/IVD
Colour	transparent
Nominal Capacity	0,2ml
Material Characteristics	polypropylene (PP)
Temperature range	-40°C - 120°C
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 6.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 125
	pcs per carton: 3.125
	pcs per pallet: 62.500
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable

Approved by nerbe plus quality department

00-000-0001/02.2020

