

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>04-072-0525</b>	<b>Revision: 9</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Deep well plate PP, premium surface, 96x2,0ml, round wells, U-shape, transparent, np pcr ready, max. RCF 4.000g, autocl. up to 121°C
<b>Dimensions</b>	ANSI/SLAS 1 & 4
	length: 127,76 mm (+/- 0,5 mm)
	width: 85,48 mm (+/- 0,5 mm)
<b>Colour</b>	transparent
<b>Nominal Capacity</b>	96x2,0ml
<b>Material Characteristics</b>	polypropylene (PP)
<b>Temperature range</b>	-80°C - 121°C
<b>Autoclavability</b>	autocl. up to 121°C
<b>Centrifugation</b>	max. RCF 4.000g
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
<b>Purity Grade</b>	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
<b>Shelf - Life</b>	not limited at controlled storage conditions
<b>Packaging</b>	pcs per bag: 5
	pcs per carton: 50
	pcs per pallet: 2.100
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C

00-000-0001/06.2020



# Technical Datasheet

Date: 02.06.22

Page 2 of 2

pal. not stackable
--------------------

<b>Approved by nerbe plus quality department</b>
--