

Technical Datasheet

		EN ISO 13485 certified
Article - No:	04-072-0565	Revision: 6
Remarks		single use only

Specification

Product	Deep well plate PP, premium surface, 96x2,0ml, square wells, V-shape, transparent, np pcr ready, max. RCF 4.000g, autocl. up to 121°C
Dimensions	ANSI/SLAS 1 & 4
	length: 127,76 mm (+/- 0,5 mm)
	width: 85,48 mm (+/- 0,5 mm)
Colour	transparent
Nominal Capacity	96x2,0ml
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 121°C
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 4.000g
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 5
	pcs per carton: 50
	pcs per pallet: 2.100
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C

Technical Datasheet

Date: 02.06.22

Page 2 of 2

pal. not stackable

Approved by nerbe plus quality department
--