

Technical Datasheet

		EN ISO 13485 certified
Article - No:	04-072-5430	Revision: 3
Remarks		single use only

Specification

Product	Reservoir PS, 25ml, high profile, 1 bottom, universal purpose, white, grad., np pcr ready, sterile R
Dimensions	length: 146,1 mm (+/- 0,2 mm)
	width: 60,5 mm (+/- 0,2 mm)
	height: 24,1 mm (+/- 0,2 mm)
Colour	white
Graduation	grad.
Nominal Capacity	25ml
Material Characteristics	polystyrene (PS)
Temperature range	-20°C - 60°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	sterile R SAL 10-6
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 5
	pcs per carton: 100
	pcs per pallet: 5.500
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C



Technical Datasheet

Date: 22.06.20

Page 2 of 2

pal. not stackable

Approved by nerbe plus quality department
--