

# Technical Datasheet

EN ISO 13485 certified		
<b>Article - No:</b>	<b>04-093-0150</b>	<b>Revision: 6</b>
<b>Remarks</b>	<b>suitable plate: 04-072-0510</b>	<b>single use only</b>

## Specification

<b>Product</b>	Silicone sealing mat, for 0,6 ml 96 deep well plates, round wells, natural, np pcr ready, autocl. up to 121°C
<b>Dimensions</b>	ANSI/SLAS 1
<b>Colour</b>	natural
<b>Material Characteristics</b>	silicone (SI)
<b>Temperature range</b>	-80°C - 121°C
<b>Autoclavability</b>	autocl. up to 121°C
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	flexible and soft
<b>Purity Grade</b>	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
<b>Shelf - Life</b>	not limited at controlled storage conditions
<b>Packaging</b>	pcs per bag: 10
	pcs per carton: 50
	pcs per pallet: 5.850
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
<b>Approved by nerbe plus quality department</b>	

00-000-0001/06.2020

