

Technical Datasheet

		EN ISO 13485 certified
Article - No:	04-095-0008	Revision: 5
Remarks		single use only

Specification

Product	Adhesive sealing film, for Realtime-PCR (qPCR), 8-strip, 50µm, -40°C up to +120°C, PE,, clear, np pcr ready
Dimensions	length: 99,00 mm (+/- 0,2 mm)
	width: 80,00 mm (+/- 0,2 mm) – 1 strip = 10,00 mm
Colour	clear
Material Characteristics	film: polyethylene (PE)
	adhesive: acrylic
Temperature range	-40°C - 120°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Sterilisation Process	N/A
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per inner cs: 400
	pcs per carton: 400
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

