

Technical Datasheet

		EN ISO 13485 certified
Article - No:	04-095-0050	Revision: 3
Remarks		single use only

Specification

Product	Adhesive sealing film, for Realtime-PCR (qPCR), 70µm, -40°C up to +104°C, PE, non-tacky, clear, np pcr ready
Dimensions	width: 77,80 mm (+/- 0,2 mm)
	length: 144,30 mm (+/- 0,2 mm)
Colour	clear
Material Characteristics	film: polyethylene (PE)
	adhesive: acrylic
Temperature range	-40°C - 104°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per inner cs: 100
	pcs per carton: 500
	pcs per pallet: 700.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable

Approved by nerbe plus quality department