

Technical Datasheet

Date: 22.06.20

Page 1 of 1

| | | |
|----------------------|--------------------|------------------------|
| | | EN ISO 13485 certified |
| Article - No: | 04-095-0200 | Revision: 8 |
| Remarks | | single use only |

Specification

| | |
|--|--|
| Product | Adhesive sealing film, for cold-storage, 60µm, -80°C up to +120°C, aluminium, pierceable, chemically resistant, np pcr ready |
| Dimensions | Length: 129mm |
| | Width 78,3mm |
| Material Characteristics | aluminium (AL) |
| Temperature range | -80°C - 120°C |
| Product Characteristics | DMSO resistant |
| Purity Grade | free of RNase and DNase |
| | free of pyrogens |
| Shelf - Life | N/A |
| Packaging | pcs per inner cs: 100 |
| | pcs per carton: 100 |
| | pcs per pallet: 700.000 |
| Storage | dry condition: 40%-60% rel. humidity |
| | bags / cases orig. closed |
| | storage temp. +2°C to max. +30°C |
| | pal. not stackable |
| Approved by nerbe plus quality department | |