

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>04-096-6008</b>	<b>Revision: 8</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Adhesive sealing film, for cold-storage, 6-strip, 50µm, -80°C up to +120°C, alu., pierceable, np pcr ready
<b>Dimensions</b>	length: 101,00 mm (+/- 0,2 mm)
	width: 60,00 mm (+/- 0,2 mm) – 1 strip = 10,00 mm
<b>Material Characteristics</b>	aluminium (AL)
<b>Temperature range</b>	-80°C - 120°C
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
<b>Purity Grade</b>	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
<b>Shelf - Life</b>	36 months at controlled storage conditions
<b>Packaging</b>	pcs per inner cs: 300
	pcs per carton: 300
	pcs per pallet: 700.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
<b>Approved by nerbe plus quality department</b>	

00-000-0001/06.2020

