

Technical Datasheet

Date: 24.05.17

Page 1 of 1

		EN ISO certified
Article - No:	04-097-5050	Revision: 4
Remarks		single use only

Specification

Product	Adhesive sealing film, permeable (tissue culture), 50µm, -4°C up to +40°C, PE, clear, np pcr ready, sterile R
Dimensions	length between perforation: 153,00 mm (+/- 0,2 mm)
	width: 82,30 mm (+/- 0,2 mm)
Colour	clear
Material Characteristics	film: polyethylene (PE)
	adhesive: acrylic
Temperature range	-4°C - 40°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
	free of endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per inner cs: 50
	pcs per carton: 500
	pcs per pallet: 3.500.000
Storage	dry condition: 40%-60% rel. humidity
	bags & cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

