

Technical Datasheet

		EN ISO 13485 certified
Article - No:	04-097-5060	Revision: 5
Remarks		single use only

Specification

Product	Adhesive sealing film, permeable (tissue culture), 140µm, -20°C up to +80°C, Viscose, natural, np pcr ready, sterile R
Dimensions	length: 143,00 mm (+/- 0,2 mm)
	width: 82,70 mm (+/- 0,2 mm)
Colour	natural
Material Characteristics	viscose (VI)
Temperature range	-20°C - 80°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	24 months at controlled storage conditions
Packaging	pcs per inner cs: 50
	pcs per carton: 50
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable

Approved by nerbe plus quality department

00-000-0001/06.2020

