

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>04-212-1203</b>	<b>Revision: 5</b>
<b>Remarks</b>	<b>piercable cap</b>	<b>single use only</b>

## Specification

<b>Product</b>	np Safelock-Cap microcentrifuge tube PP, 1,5ml, attached cap, boil-proof, yellow, grad., writ. area, np pcr ready, max. RCF 30.000g, autocl. up to 121°C
<b>Dimensions</b>	outer diameter: 10,80 mm (+/- 0,2 mm)
	height (without cap): 38,90 mm (+/- 0,2 mm)
<b>Colour</b>	yellow
<b>Graduation</b>	grad., writ. area
<b>Nominal Capacity</b>	1,5ml
<b>Material Characteristics</b>	polypropylene (PP)
<b>Temperature range</b>	-70°C - 75°C
	in waterbath up to 120 °C
<b>Autoclavability</b>	autocl. up to 121°C
<b>Centrifugation</b>	max. RCF 30.000g
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	constantly usable at temperature of 80 °C up to 100 °C
<b>Purity Grade</b>	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
<b>Shelf - Life</b>	not limited at controlled storage conditions
<b>Packaging</b>	pcs per bag: 1.000
	pcs per carton: 8.000
	pcs per pallet: 160.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable

# Technical Datasheet

Date: 03.06.22

Page 2 of 2



EN ISO 13485 certified

Approved by nerbe plus quality department

nerbe plus GmbH & Co. KG  
Porschestr. 25  
D-21423 Winsen/Luhe  
Tel: +49 (0)4171 7849 0  
Fax: +49 (0)4171 6 45 45

info@nerbe-plus.de  
www.nerbe-plus.de  
HRA Lüneburg 203240  
USt.-ID: DE 32583 0652

Pers. haftende Gesellschafterin:  
NERBE R&D GmbH, Winsen/Luhe  
GF: Philipp J.F. Nerbe  
HRB Lüneburg 208826

Deutsche Bank AG  
Kto. 0204 41200  
IBAN: DE41 2407 0075 0020 4412 00  
SWIFT / BIC DEUTDE2H240

Hamburger Sparkasse  
Kto. 1387 140831  
IBAN: DE56 2005 0550 1387 1408 31  
SWIFT / BIC HASPDE33XXX

