

Technical Datasheet

		EN ISO 13485 certified
Article - No:	04-232-3100	Revision: 6
Remarks	piercable cap	single use only

Specification

Product	Microcentrifuge tube PP, 2ml, attached cap, transparent, grad., writ. area, np pcr ready, max. RCF 18.000g, autocl. up to 121°C
Colour	transparent
Graduation	grad., writ. area
Nominal Capacity	2ml
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 121°C
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 18.000g
Product Characteristics	biologically inert
	physiologically harmless
	not heat resistant
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Packaging	pcs per bag: 200
	pcs per carton: 5.000
	pcs per pallet: 140.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

