

Technical Datasheet

		EN ISO 13485 certified
Article - No:	04-332-3000	Revision: 6
Remarks	piercable cap	single use only

Specification

Product	Microcentrifuge tube PP, 0,5ml, attached cap, heat-resistant, transparent, grad., writ. area, np pcr ready, max. RCF 20.000g, autocl. up to 121°C
Colour	transparent
Graduation	grad., writ. area
Nominal Capacity	0,5ml
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 121°C
Autoclavability	autocl. up to 121°C
Centrifugation	max. RCF 20.000g
Product Characteristics	biologically inert
	physiologically harmless
	constantly usable at temperature of 80 °C up to 100 °C
Purity Grade	free of PCR inhibitors
	free of human DNA und RNA
	free of DNase and RNase
	free of pyrogens
	free of endotoxins
Packaging	pcs per bag: 200
	pcs per carton: 5.000
	pcs per pallet: 240.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +2°C to max. +30°C
	pal. not stackable

Approved by nerbe plus quality department

00-000-0001/06.2020

