

Technical Datasheet

		EN ISO 13485 certified
Article - No:	05-041-0000	Revision: 2
Remarks		single use only

Specification

Product	Analyzer cup PS, 0,5ml, Ø13,55x24,5mm, for RA 1000 from Technikon, Gemsac, Kone and Astra, transparent, CE/IVD
Classification	CE/IVD
Dimensions	diameter: 13,55 mm (+/- 0,2 mm) height: 24,50 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	0,5ml
Material Characteristics	polystyrene (PS)
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62 Food contact: raw material in compliance with regulation EC 1935/2004
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 1.000 pcs per carton: 10.000 pcs per pallet: 120.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	