

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>05-500-1015</b>	<b>Revision: 3</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Micro screw cap tube PP, 1,5ml, assembeled screw cap PP with o-ring, tube natural, cap transparent sterile R, max. RCF 20.000g
<b>Colour</b>	tube natural, cap transparent
<b>Nominal Capacity</b>	1,5ml
<b>Material Characteristics</b>	polypropylene (PP)
<b>Temperature range</b>	-196°C - 121°C
<b>Centrifugation</b>	max. RCF 20.000g
<b>Purity Grade</b>	free of RNase and DNase free of pyrogens
<b>Sterilisation Process</b>	radiation sterilization SAL 10-6
<b>Packaging</b>	pcs per bag: 100 pcs per inner cs: 500 pcs per carton: 4.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable
<b>Approved by nerbe plus quality department</b>	

00-000-0001/06.2020

