

Technical Datasheet

Date: 09.06.22

Page 1 of 1

		EN ISO 13485 certified
Article - No:	05-501-0020	Revision: 3
Remarks		single use only

Specification

Product	Micro screw cap tube PP, skirted, 2ml, with enclosed screw cap PP with O-ring, tube natural, cap transparent, max. RCF 20.000g
Colour	tube natural, cap transparent
Nominal Capacity	2ml
Material Characteristics	polypropylene (PP)
Temperature range	-196°C - 121°C
Centrifugation	max. RCF 20.000g
Purity Grade	free of RNase and DNase free of pyrogens
Packaging	pcs per bag: 500 pcs per carton: 4.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	