

Technical Datasheet

		EN ISO 13485 certified
Article - No:	05-501-1015	Revision: 4
Remarks		single use only

Specification

Product	Micro screw cap tube PP, skirted, 1,5ml, assembeled screw cap PP with o-ring, tube natural, cap transparent, sterile R, max. RCF 20.000g
Colour	tube natural, cap transparent
Nominal Capacity	1,5ml
Material Characteristics	polypropylene (PP)
Temperature range	-196°C - 121°C
Centrifugation	max. RCF 20.000g
Purity Grade	free of RNase and DNase free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Packaging	pcs per bag: 100 pcs per carton: 4.000 pcs per pallet: 96.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

