

Technical Datasheet

EN ISO 13485 certified	
Article - No:	05-733-1000
Remarks	Revision: 2 single use only

Specification

Product	Tubes PE for scintillation bottles, 4ml, enclosed screw cap PP, Ø16x57 mm, natural
Dimensions	diameter: 16,00 mm (+/- 0,2 mm)
	height: 57,00 mm (+/- 0,2 mm)
Colour	natural
Nominal Capacity	4ml
Material Characteristics	polyethylene (PE)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 1.000
	pcs per carton: 2.000
	pcs per pallet: 84.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

