

Technical Datasheet

		EN ISO 13485 certified
Article - No:	05-900-1005	Revision: 2
Remarks		single use only

Specification

Product	Cryo vial PP, screw cap PE, 0,5ml, self-standing, transparent, grad., sterile R, autocl.
Dimensions	height: 28.6mm
Colour	transparent
Graduation	grad.
Nominal Capacity	0,5ml
Material Characteristics	tube: polypropylene (PP) cap: polyethylene (PE)
Temperature range	-196°C - 121°C
Autoclavability	autocl.
Purity Grade	free of PCR inhibitors free of human DNA free of RNase and DNase free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Packaging	pcs per bag: 100 pcs per carton: 1.000 pcs per pallet: 99.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	