

Technical Datasheet

		EN ISO 13485 certified
Article - No:	05-900-1010	Revision: 3
Remarks		single use only

Specification

Product	Cryo vial PP, screw cap PE, self-standing, transparent, grad., sterile R, autocl.
Dimensions	height: 36,1mm
Colour	transparent
Graduation	grad.
Material Characteristics	tube: polypropylene (PP) cap: polyethylene (PE)
Temperature range	-196°C - 121°C
Autoclavability	autocl.
Purity Grade	free of PCR inhibitors free of human DNA free of RNase and DNase free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Packaging	pcs per bag: 100 pcs per carton: 1.000 pcs per pallet: 99.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

