

# Technical Datasheet

Date: 28.06.22

Page 1 of 1

		EN ISO 13485 certified
<b>Article - No:</b>	<b>05-902-0007</b>	<b>Revision: 3</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Insert for cryo tube lid, white, autocl.
<b>Colour</b>	white
<b>Temperature range</b>	-196°C - 121°C
<b>Autoclavability</b>	autocl.
<b>Packaging</b>	pcs per bag: 500 pcs per carton: 500
<b>Storage</b>	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable
<b>Approved by nerbe plus quality department</b>	