

Technical Datasheet

Date: 28.06.22

Page 1 of 1

		EN ISO 13485 certified
Article - No:	05-902-0009	Revision: 5
Remarks		single use only

Specification

Product	Insert for cryo tube lid, 5 different colors, autocl.
Colour	5 different colors
Temperature range	-196°C - 121°C
Autoclavability	autocl.
Packaging	pcs per bag: 500 pcs per carton: 500
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +2°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	