

Technical Datasheet

		EN ISO 13485 certified
Article - No:	06-361-2008	Revision: 14
Remarks		single use only

Specification

Product	Pipette tips PP, premium surface, 0,1-20µl, extra long, transparent, moulded rings, np premium bio, sterile R, autocl. up to 121°C
Dimensions	length: 45,75 mm (+/- 0,2 mm)
	inner diameter: 0,40 mm (+/- 0,2 mm)
Colour	transparent
Graduation	moulded rings
Nominal Capacity	0,1-20µl
Material Characteristics	polypropylene (PP), premium surface
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	sterile R SAL 10-6
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of ATP
	free of pyrogens and endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 1.000
	pcs per carton: 20.000
	pcs per pallet: 840.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable



Technical Datasheet

Date: 02.06.22

Page 2 of 2

Approved by nerbe plus quality department