

Technical Datasheet

		EN ISO 13485 certified
Article - No:	06-602-5300	Revision: 8
Remarks		single use only

Specification

Product	Filter tip PP, premium surface, 0,1-10µl, short, transparent, moulded rings, np premium bio, sterile R
Dimensions	length: 31,30 mm (+/- 0,2 mm)
	inner diameter: 0,50 mm (+/- 0,2 mm)
Colour	transparent
Graduation	moulded rings
Nominal Capacity	0,1-10µl
Material Characteristics	polypropylene (PP), premium surface
	filter: high density polyethylene (PE)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	sterile R SAL 10-6
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of ATP
	free of pyrogens and endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 960
	pcs per carton: 5.760
	pcs per pallet: 172.800
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable



Technical Datasheet

Date: 02.06.22

Page 2 of 2

Approved by nerbe plus quality department