

Technical Datasheet

		EN ISO 13485 certified
Article - No:	06-672-5300	Revision: 7
Remarks		single use only

Specification

Product	Filter tip PP, premium surface, 0-300µl, 96 pcs/rack, natural, moulded rings, np premium bio, sterile R
Dimensions	length: 59,35 mm (+/- 0,2 mm)
	inner diameter: 0,52 mm (+/- 0,2 mm)
Colour	natural
Graduation	moulded rings
Nominal Capacity	0-300µl
Material Characteristics	polypropylene (PP), premium surface
	filter: high density polyethylene (PE)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	sterile R SAL 10-6
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of ATP
	free of pyrogens and endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 960
	pcs per carton: 5.760
	pcs per pallet: 172.800
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable

00-000-0001/06.2020



Technical Datasheet

Date: 02.06.22

Page 2 of 2

Approved by nerbe plus quality department