

Technical Datasheet

Date: 02.06.22

Page 1 of 1

EN ISO 13485 certified		
Article - No:	07-012-8100	Revision: 2
Remarks	Suitable for: Multipette® 4780, ripette®, HandyStep®, Minilab 100/101 STEPMATE, EasyStep, Distriman®, HandyStep®electronic	single use only

Specification

Product	Dispenser-tips PP/PE, 0,5ml, single peel-pack, grad., np pcr ready, sterile
Dimensions	length: 100,50 mm (+/- 0,2 mm)
Graduation	grad.
Nominal Capacity	0,5ml
Material Characteristics	cylinder: polypropylene (PP) piston: polyethylene (PE)
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors free of human DNA und RNA free of DNase and RNase free of pyrogens free of endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per carton: 1.200 pcs per pallet: 9.600 pcs per dsp case: 100
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	