

# Technical Datasheet

EN ISO 13485 certified		
<b>Article - No:</b>	<b>07-032-0100</b>	<b>Revision: 3</b>
<b>Remarks</b>	<b>Suitable for: Multipette® 4780, ripette®, HandyStep®, Minilab 100/101 STEPMATE, EasyStep, Distriman®, HandyStep®electronic</b>	<b>single use only</b>

## Specification

<b>Product</b>	Dispenser-tips PP/PE, 2,5ml, grad., np pcr ready
<b>Dimensions</b>	length: 107,00 mm (+/- 0,2 mm)
<b>Graduation</b>	grad.
<b>Nominal Capacity</b>	2,5ml
<b>Material Characteristics</b>	cylinder: polypropylene (PP) piston: polyethylene (PE)
<b>Product Characteristics</b>	biologically inert physiologically harmless no heavy metals according to EC/94/62
<b>Purity Grade</b>	free of PCR inhibitors free of human DNA und RNA free of DNase and RNase free of pyrogens free of endotoxins
<b>Shelf - Life</b>	60 months at controlled storage conditions
<b>Packaging</b>	pcs per carton: 4.800 pcs per pallet: 38.400 pcs per dsp case: 200
<b>Storage</b>	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
<b>Approved by nerbe plus quality department</b>	