

Technical Datasheet

EN ISO 13485 certified		
Article - No:	07-062-0100	Revision: 3
Remarks	Suitable for: Multipette® 4780, ripette®, HandyStep®, Minilab 100/101 STEPMATE, EasyStep, Distriman®, HandyStep®electronic	single use only

Specification

Product	Dispenser-tips PP/PE, 25ml, grad., np pcr ready
Dimensions	length: 126,00 mm (+/- 0,2 mm)
Graduation	grad.
Nominal Capacity	25ml
Material Characteristics	cylinder: polypropylene (PP) piston: polyethylene (PE)
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors free of human DNA und RNA free of DNase and RNase free of pyrogens free of endotoxins
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per carton: 800 pcs per pallet: 6.400 pcs per dsp case: 50
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	