

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>07-112-0000</b>	<b>Revision: 3</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Pipette tip PP "classic", 0,1-20µl, universal fit (Eppendorf, Biohit, Socorex), transparent, np pcr ready
<b>Dimensions</b>	length: 45,50 mm (+/- 0,2 mm)
<b>Colour</b>	transparent
<b>Nominal Capacity</b>	0,1-20µl
<b>Material Characteristics</b>	polypropylene (PP)
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
<b>Purity Grade</b>	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
<b>Shelf - Life</b>	60 months at controlled storage conditions
<b>Packaging</b>	pcs per bag: 1.000
	pcs per carton: 8.000
	pcs per pallet: 864.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
<b>Approved by nerbe plus quality department</b>	

00-000-0001/06.2020

