

Technical Datasheet

EN ISO 13485 certified		
Article - No:	07-122-0073	Revision: 2
Remarks	Suitable for common pipettes e.g. Eppendorf, Biohit, Gilson or Brand	single use only

Specification

Product	Pipette tip PP "classic", 1-200µl, universal fit (Eppendorf, Biohit, Socorex), yellow, np pcr ready, autocl. up to 121°C
Dimensions	length: 52,60 mm (+/- 0,2 mm)
Colour	yellow
Nominal Capacity	1-200µl
Material Characteristics	polypropylene (PP)
Temperature range	-20°C - +90°C
Autoclavability	autocl. up to 121°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Packaging	pcs per bag: 1.000
	pcs per carton: 25.000
	pcs per pallet: 500.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	