

Technical Datasheet

EN ISO 13485 certified		
Article - No:	07-132-0095	Revision: 6
Remarks	Suitable for common pipettes e.g. Eppendorf, Biohit, Gilson or Brand	single use only

Specification

Product	Pipette tip PP "classic", 100-1.000µl, universal fit (Eppendorf, Biohit, Socorex), blue, autocl. up to 121°C
Dimensions	inner diameter conus: 7,70 mm (+/- 0,2 mm)
	length: 70,40 mm (+/- 0,2 mm)
Colour	blue
Nominal Capacity	100-1.000µl
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 121°C
Autoclavability	autocl. up to 121°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Packaging	pcs per bag: 1.000
	pcs per carton: 10.000
	pcs per pallet: 200.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

