

Technical Datasheet

		EN ISO 13485 certified
Article - No:	07-360-2015	Revision: 6
Remarks		single use only

Specification

Product	Pipette tips PP, premium surface, 0,1-10µl, short transparent, moulded rings, np pcr ready, autocl. up to 121°C
Dimensions	length: 31,20 mm (+/- 0,2 mm) inner diameter: 0,46 mm (+/- 0,2 mm)
Colour	transparent
Graduation	moulded rings
Nominal Capacity	0,1-10µl
Material Characteristics	polypropylene (PP), premium surface
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62 sterile SAL 10-3
Purity Grade	free of RNase and Dnase free of pyrogens free of PCR inhibitors free of human DNA free of endotoxins free of ATP
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 960 pcs per carton: 4.800 pcs per pallet: 201.600
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable

Technical Datasheet

Date: 03.06.22

Page 2 of 2

Approved by nerbe plus quality department