

Technical Datasheet

		EN ISO 13485 certified
Article - No:	07-366-2015	Revision: 5
Remarks		single use only

Specification

Product	Pipette tips PP, premium surface, 0,5-200µl, extra long, 96 pcs/rack, transparent, moulded rings, np pcr ready, autocl. up to 121°C
Dimensions	length: 59,54 mm (+/- 0,2 mm)
	inner diameter: 0,51 mm (+/- 0,2 mm)
Colour	transparent
Graduation	moulded rings
Nominal Capacity	0,5-200µl
Material Characteristics	polypropylene (PP), premium surface
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	sterile SAL 10-3
Purity Grade	free of RNase and Dnase
	free of pyrogens
	free of PCR inhibitors
	free of human DNA
	free of endotoxins
	free of ATP
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 960
	pcs per carton: 4.800
	pcs per pallet: 201.600
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable



Technical Datasheet

Date: 03.06.22

Page 2 of 2

Approved by nerbe plus quality department